

Comparison between MiSeq and HiSeq using Paired Data

Quanhui Sheng, Bojana Jovanovic, Scott Austin Beeler

March 4, 2014

Contents

| | | |
|----------|--|-----------|
| 1 | Purpose | 2 |
| 2 | Material and method | 2 |
| 2.1 | Dataset | 2 |
| 2.2 | Raw file quality control | 2 |
| 2.3 | Terminal 'N' trimming | 2 |
| 2.4 | Mapping | 2 |
| 2.5 | Counting genes | 2 |
| 2.6 | FPKM of genes | 3 |
| 3 | Result | 3 |
| 3.1 | Raw file quality control | 3 |
| 3.2 | Mapping quality | 10 |
| 3.3 | Correlation of count between MiSeq and HiSeq | 10 |
| 3.4 | Correlation of FPKM between MiSeq and HiSeq | 10 |
| 3.5 | Heatmap | 11 |
| 3.6 | PCA | 12 |
| 3.7 | Differential expression comparison | 13 |
| 4 | Discussion | 13 |

1 Purpose

The MiSeq and HiSeq sequencing machines have been widely used in next generation sequencing projects. Here, we want to validate the possibility to combine the data from those two machines together for further analysis.

2 Material and method

2.1 Dataset

Three samples were pair-end sequenced by both MiSeq and HiSeq machines, as table 1 illustrated.

Table 1: Sample Information

| Type | Source | Alias |
|-------|-----------------|--------------|
| HiSeq | 2059-JP-6 | HiSeqSample1 |
| HiSeq | 2059-JP-7 | HiSeqSample2 |
| HiSeq | 2059-JP-9 | HiSeqSample3 |
| MiSeq | IG-33_408637_S4 | MiSeqSample1 |
| MiSeq | IG-34_S3 | MiSeqSample2 |
| MiSeq | IG-39_408648_S5 | MiSeqSample3 |

2.2 Raw file quality control

FastQC v0.10.1 was used to evaluate the quality of FastQ files.

```
fastqc -t 2 -o HiSeqSample1 2059-JP-6_1.fastq.gz 2059-JP-6_2.fastq.gz
```

2.3 Terminal 'N' trimming

The terminal 'N' in data was trimmed by in-house software cqstools.

```
mono-sgen CQS.Tools.exe fastq.trimmer -n -z -i 2059-JP-6_1.fastq.gz -o 2059-JP-6_1_trim.fastq.gz
```

2.4 Mapping

TopHat v2.0.9 was used to map reads to human genome 19.

```
tophat2 --segment-length 25 -r 0 -p 8 --keep-fastq-order --no-coverage-search \
--rg-id HiSeqSample1 --rg-sample HiSeqSample1 --rg-library HiSeqSample1 \
-G Homo_sapiens.GRCh37.73.gtf --transcriptome-index=hg19.GRCh37.73 \
-o . bowtie2_index/hg19 2059-JP-6_1.fastq.gz 2059-JP-6_2.fastq.gz
```

```
mv accepted_hits.bam HiSeqSample1.bam
```

```
mv accepted_hits.bam.bai HiSeqSample1.bam.bai
```

2.5 Counting genes

HTSeq v0.5.4p3 was used to count genes.

```
samtools sort -n -@ 8 HiSeqSample1.bam HiSeqSample1.sortedname
```

```
samtools view HiSeqSample1.sortedname.bam | htseq-count -q -m intersection-nonempty \
-s no -i gene_id - Homo_sapiens.GRCh37.73.gtf >HiSeqSample1.count
```

2.6 FPKM of genes

Cuffdiff v2.1.1 were used to estimate gene FPKM value. cqstools was used to extract FPKM values.









```
cuffdiff -p 8 -u -N -o . -L MiSeq,HiSeq -b hg19_chr.fa Homo.sapiens.GRCh37.73.gtf \  
  MiSeqSample1.bam,MiSeqSample2.bam,MiSeqSample3.bam \  
  HiSeqSample1.bam,HiSeqSample2.bam,HiSeqSample3.bam  
  
cqstools cuffdiff_count -i genes.read_group_tracking -s gene_exp.diff \  
  -m 20140218_bojana_MiSeq_HiSeq_group_sample.map -o untrim_cuffdiff
```









3 Result









3.1 Raw file quality control









No adapter was found. There are three questionable observations:









1. The GC contribution of HiSeq data was shifted to 0.4.
2. There is a 'N' in 3' of each HiSeq reads.
3. There are a lot of reads in MiSeq have no GC, which was caused by the reads whose sequence was all 'N'.





















3.2 Mapping quality

We used cqstools to trim the terminal 'N' of the reads. Both pretrim reads (table 2) and posttrim reads (table 3) were mapped to genome and summerized.

Table 2: Mapping summary of untrimmed reads

| | LeftReads | LMapped | LRate | RightReads | RMapped | RRate | AlignedPairs | ADiscordant |
|--------|-----------|----------|-------|------------|----------|-------|--------------|-------------|
| HiSeq1 | 24712107 | 21215897 | 85.9% | 24712107 | 19631934 | 79.4% | 19103150 | 6.3% |
| HiSeq2 | 24420179 | 19812406 | 81.1% | 24420179 | 17402588 | 71.3% | 16974487 | 6.3% |
| HiSeq3 | 25989600 | 22146933 | 85.2% | 25989600 | 19273082 | 74.2% | 18718125 | 6.4% |
| MiSeq1 | 6270058 | 4915197 | 78.4% | 6270058 | 4746878 | 75.7% | 4444182 | 6.0% |
| MiSeq2 | 7402812 | 6548949 | 88.5% | 7402812 | 6414037 | 86.6% | 6307380 | 2.5% |
| MiSeq3 | 6338368 | 5162507 | 81.4% | 6338368 | 4872895 | 76.9% | 4604894 | 4.2% |

Table 3: Mapping summary of trimmed reads

| | LeftReads | LMapped | LRate | RightReads | RMapped | RRate | AlignedPairs | ADiscordant |
|--------|-----------|----------|-------|------------|----------|-------|--------------|-------------|
| HiSeq1 | 24712107 | 22637599 | 91.6% | 24712107 | 21039704 | 85.1% | 20492712 | 8.5% |
| HiSeq2 | 24420179 | 21253887 | 87.0% | 24420179 | 18784680 | 76.9% | 18344411 | 8.7% |
| HiSeq3 | 25989600 | 23550896 | 90.6% | 25989600 | 20593996 | 79.2% | 19994182 | 7.8% |
| MiSeq1 | 5854210 | 4909295 | 83.9% | 5856888 | 4731464 | 80.8% | 4101953 | 98.2% |
| MiSeq2 | 7345681 | 6544854 | 89.1% | 7346956 | 6405153 | 87.2% | 5787518 | 98.8% |
| MiSeq3 | 6045508 | 5157791 | 85.3% | 6046582 | 4859000 | 80.4% | 4271930 | 98.9% |

Trimming terminal 'N' increased the mapping sensitivity of HiSeq data but introduced a lot of discordant mapped pairs in MiSeq data. So, we chose the posttrim HiSeq and pretrim MiSeq mapping result in the following analysis.

3.3 Correlation of count between MiSeq and HiSeq

Spearman correlation was calculated between MiSeq and HiSeq gene count for each sample (table 4).

Table 4: Spearman Correlation of Count between MiSeq and HiSeq

| sample | spcorr |
|---------|--------|
| Sample1 | 0.71 |
| Sample2 | 0.74 |
| Sample3 | 0.86 |

3.4 Correlation of FPKM between MiSeq and HiSeq

Spearman correlation was calculated between MiSeq and HiSeq gene FPKM value for each sample (table 5).


Table 5: Spearman Correlation of FPKM between MiSeq and HiSeq

| sample | spcorr |
|---------|--------|
| Sample1 | 0.62 |
| Sample2 | 0.64 |
| Sample3 | 0.77 |

3.5 Heatmap


Variance stabilizing transformed (VST) value of each gene was calculated by DESeq2 from gene count. The 500 genes with largest interquartile range (IQR) of VST value were selected for drawing heatmap.

Based on heatmap, the MiSeq and HiSeq data from sample 3 were clustered together as we expected. The spearman correlation coefficient of sample 3 (0.86) between MiSeq and HiSeq was also higher than sample 1 and 2 (0.71 and 0.74).



3.6 PCA

The VST values of all genes were used in PCA analysis. Based on PCA image, the MiSeq and HiSeq data from sample 3 were also more similar than the data from sample 1 and 2.



3.7 Differential expression comparison

Ideally, there should be no difference between MiSeq and HiSeq paired data. But actually, DESeq2 detected 6787 genes differential expressed with adjust pvalue less than 0.05.

Table 6: Top 50 differential expressed genes

| Name | M1 | M2 | M3 | H1 | H2 | H3 | log2(fc) | pvalue | padj |
|----------|------|------|------|-----|-----|------|----------|--------|------|
| SEPT9 | 602 | 831 | 612 | 38 | 83 | 94 | -4.42 | 0.00 | 0.00 |
| COL6A2 | 944 | 2257 | 646 | 151 | 425 | 216 | -3.34 | 0.00 | 0.00 |
| SREBF1 | 445 | 1423 | 1640 | 60 | 181 | 388 | -3.77 | 0.00 | 0.00 |
| ATN1 | 1427 | 792 | 626 | 123 | 100 | 149 | -3.97 | 0.00 | 0.00 |
| DNAJB2 | 216 | 318 | 348 | 12 | 23 | 43 | -4.70 | 0.00 | 0.00 |
| SNRNP70 | 775 | 858 | 795 | 92 | 145 | 160 | -3.75 | 0.00 | 0.00 |
| MARK2 | 578 | 754 | 492 | 116 | 159 | 190 | -3.11 | 0.00 | 0.00 |
| ATXN2L | 391 | 543 | 538 | 60 | 89 | 181 | -3.41 | 0.00 | 0.00 |
| ACTN4 | 2335 | 2482 | 2192 | 626 | 695 | 1047 | -2.74 | 0.00 | 0.00 |
| PLEC | 1211 | 1587 | 365 | 98 | 221 | 52 | -4.15 | 0.00 | 0.00 |
| SCRIB | 454 | 524 | 251 | 27 | 61 | 37 | -4.35 | 0.00 | 0.00 |
| NBEAL2 | 221 | 556 | 314 | 25 | 44 | 50 | -4.22 | 0.00 | 0.00 |
| COL18A1 | 148 | 615 | 96 | 15 | 61 | 19 | -4.07 | 0.00 | 0.00 |
| ADRBK1 | 340 | 273 | 473 | 52 | 40 | 150 | -3.47 | 0.00 | 0.00 |
| FASN | 835 | 424 | 2844 | 117 | 43 | 730 | -3.76 | 0.00 | 0.00 |
| PTK7 | 425 | 820 | 736 | 86 | 157 | 301 | -3.11 | 0.00 | 0.00 |
| KDM5C | 545 | 499 | 867 | 87 | 123 | 280 | -3.21 | 0.00 | 0.00 |
| PHRF1 | 154 | 228 | 337 | 16 | 17 | 49 | -4.31 | 0.00 | 0.00 |
| SBF1 | 206 | 337 | 158 | 22 | 56 | 33 | -3.75 | 0.00 | 0.00 |
| AP3D1 | 315 | 1038 | 487 | 59 | 188 | 205 | -3.16 | 0.00 | 0.00 |
| SF3A2 | 190 | 702 | 213 | 28 | 76 | 40 | -3.88 | 0.00 | 0.00 |
| EPN1 | 171 | 212 | 259 | 22 | 27 | 44 | -3.90 | 0.00 | 0.00 |
| LRP5 | 179 | 198 | 230 | 13 | 17 | 49 | -4.19 | 0.00 | 0.00 |
| PRR12 | 166 | 92 | 175 | 9 | 4 | 12 | -5.10 | 0.00 | 0.00 |
| FURIN | 212 | 358 | 125 | 14 | 38 | 29 | -4.11 | 0.00 | 0.00 |
| TAOK2 | 283 | 183 | 302 | 12 | 7 | 44 | -4.88 | 0.00 | 0.00 |
| CEP170B | 321 | 223 | 473 | 21 | 18 | 111 | -4.20 | 0.00 | 0.00 |
| LMNB2 | 409 | 379 | 500 | 27 | 19 | 99 | -4.48 | 0.00 | 0.00 |
| KRT5 | 1202 | 445 | 191 | 266 | 138 | 100 | -2.73 | 0.00 | 0.00 |
| AKT1 | 515 | 768 | 3426 | 80 | 211 | 1365 | -3.06 | 0.00 | 0.00 |
| SLC2A4RG | 97 | 298 | 211 | 17 | 61 | 76 | -3.19 | 0.00 | 0.00 |
| FBRSL1 | 139 | 244 | 272 | 19 | 23 | 43 | -4.03 | 0.00 | 0.00 |
| MAPK8IP3 | 297 | 283 | 158 | 31 | 20 | 15 | -4.44 | 0.00 | 0.00 |
| ZBTB7A | 302 | 376 | 201 | 76 | 91 | 81 | -2.90 | 0.00 | 0.00 |
| EHBP1L1 | 338 | 177 | 152 | 6 | 9 | 4 | -5.78 | 0.00 | 0.00 |
| USF2 | 419 | 327 | 376 | 90 | 60 | 104 | -3.27 | 0.00 | 0.00 |
| SEMA3F | 86 | 193 | 224 | 6 | 13 | 42 | -4.33 | 0.00 | 0.00 |
| NFIX | 1051 | 829 | 50 | 182 | 214 | 21 | -3.03 | 0.00 | 0.00 |
| PPP1R14B | 269 | 669 | 347 | 33 | 150 | 121 | -3.30 | 0.00 | 0.00 |
| PRKCSH | 805 | 1117 | 791 | 205 | 232 | 375 | -2.89 | 0.00 | 0.00 |

4 Discussion

Due to unknown reason, there are huge difference between the sequencing data from MiSeq and HiSeq platforms. Combining those data together to do analysis may introduce more false positives.